


Human RAET1E / Retinoic acid early transcript 1E / UL16-binding protein 4 / NKG2DL4 / ULBP4



Product Name	Human RAET1E / Retinoic acid early transcript 1E / UL16-binding protein 4 / NKG2DL4 / ULBP4
Catalog No	IT-300-148p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human RAET1E / Retinoic acid early transcript 1E / UL16-binding protein 4 / NKG2DL4 / ULBP4 (aa 31-225; GenBank# Q8TD07) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥85%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: Human NKG2DL4 / ULBP4

Note: