


Human ULBP3 / UL16-binding protein 3 / NKG2DL3 / RAET1N



Product Name	Human ULBP3 / UL16-binding protein 3 / NKG2DL3 / RAET1N
Catalog No	IT-300-151p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human ULBP3 / UL16-binding protein 3 / NKG2DL3 / RAET1N (aa 30-217; GenBank# Q9BZM4) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: Human NKG2DL3 / RAET1N

Note: