


Human Interleukin-4 receptor subunit alpha / IL4RA / CD124 / IL4-BP



Product Name	Human Interleukin-4 receptor subunit alpha / IL4RA / CD124 / IL4-BP
Catalog No	IT-300-159p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Interleukin-4 receptor subunit alpha / IL4RA / CD124 / IL4-BP (aa 26-232; GenBank# P24394) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Interleukin-4 receptor subunit alpha / IL4RA / CD124 / IL4-BP

Note: