


Human Kallikrein-8 / KLK8 / TADG-14



Product Name	Human Kallikrein-8 / KLK8 / TADG-14
Catalog No	IT-400-024p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Kallikrein-8 / KLK8 / TADG-14 (aa 29-260; GenBank# O60259) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Kallikrein-8 / KLK8 / TADG-14

Note: