


gp140(CNE55)(HIV/Clade CRF01)



Product Name	gp140(CNE55)(HIV/Clade CRF01)
Catalog No	IT-001-RP12ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp140(CNE55)(HIV/Clade CRF01) (aa 30-659; GenBank# ADI62587) with a 6xHis tag at C-terminus. Most are in trimer form.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: gp140(CNE55)(HIV/Clade CRF01)

Note: