

# Human Mucin-16/MUC16/CA125/CA-125



Product Name	Human Mucin-16 / MUC16 / CA125 / CA-125
Catalog No	IT-400-027p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Mucin-16 / MUC16 / CA125 / CA-125 (aa 12518-13168; GenBank# Q8WXI7) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Mucin-16 / MUC16 / CA125 / CA-125

Note: