


3C-like proteinase (CoVID-19/Wuhan), His-tagged



Product Name	3C-like proteinase (CoVID-19/Wuhan), His-tagged
Catalog No	IT-002-037Ep
Source	Recombinant protein purified from E.coli culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in 50mM Tris-HCl and 150mM NaCl.
Protein	The new coronavirus 3C-like proteinase (CoVID-19/Wuhan) (aa 3264-3569; GenBank# QIZ13716) with 6xHis tag at N-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: 3C-like proteinase (CoVID-19/Wuhan), his-tagged

Note: