


Human Phosphatidylcholine-sterol acyltransferase / LCAT



Product Name	Human Phosphatidylcholine-sterol acyltransferase / LCAT
Catalog No	IT-400-028p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Phosphatidylcholine-sterol acyltransferase / LCAT (aa 25-440; GenBank# P04180) with 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Phosphatidylcholine-sterol acyltransferase / LCAT

Note: