


Human ANA/BTG3/TOB5



| | |
|----------------|--|
| Product Name | Human ANA/BTG3/TOB5 |
| Catalog No | IT-010-027Ep |
| Source | Recombinant protein purified from E.coli culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in 50mM Tris-HCl and 150mM NaCl. |
| Protein | Human ANA/BTG3/TOB5 (aa 1-252; GenBank# Q14201) with 6xHis tag at C-terminus. |
| Applications | WB standard, antibody ELISA, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

Reference



E: Human ANA/BTG3/TOB5

Note: