

Anti-Spike & NP (SARS-CoV-2/Wuhan), Rabbit Polyclonal



Product Name	Anti-Spike & NP (SARS-CoV-2/Wuhan)
Catalog No	IT-002-030SN
Source	Rabbit Polyclonal
Isotype	IgG
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing Spike (SARS-CoV-2/Wuhan)(IT-002-032p) and NP (IT-002-033Ep).
Applications	WB, ELISA, IP, IF etc.
Specificity	React to Spike and NP (SARS-CoV-2/Wuhan) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

