

Human Natural killer cells antigen CD94 / KLRD1



Product Name	Human Natural killer cells antigen CD94 / KLRD1
Catalog No	IT-300-177p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	N-terminal 10xHis tagged Human Natural killer cells antigen CD94 / KLRD1 (aa 34-179; GenBank# Q13241).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Natural killer cells antigen CD94 / KLRD1

Note: