

HA(Δ TM)(A/swine/Jiangsu/J004/2018) (H1N1)



Product Name	HA(Δ TM)(A/swine/Jiangsu/J004/2018)(H1N1)
Catalog No	IT-003-SW17 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/swine/Jiangsu/J004/2018) (H1N1)(aa 18-529; GenBank# QEM24690), a G4 EA strain.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazaIT-300-177prdous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

