

Human Low-density lipoprotein receptor class A domain-containing protein 3/ LDLRAD3



| | |
|----------------|--|
| Product Name | Human Low-density lipoprotein receptor class A domain-containing protein 3 / LDLRAD3 |
| Catalog No | IT-600-013p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | 6xHis tagged Human LDLRAD3 (aa 18-173; GenBank# Q86YD5). LDLRAD3 is the receptor for VEEV. |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

