

gp120(2759058)(HIV/Clade C)



| | |
|----------------|--|
| Product Name | gp120(2759058)(HIV/Clade C) |
| Catalog No | IT-001-RC16p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | gp120(2759058)(HIV/Clade C) (aa 33-497; GenBank# HQ615959) with a 6xHis tag at C-terminus. |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |
| Reference | |



E: gp120(2759058)(HIV/Clade C)

Note: