


## Small delta antigen/HDVAg (Hepatitis delta virus /ISOLATE AMERICAN)



Product Name	Small delta antigen (Hepatitis delta virus/ISOLATE AMERICAN)
Catalog No	IT-004-024Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged Small delta antigen (Hepatitis delta virus/ISOLATE AMERICAN) (a.a. 1-195, GenBank# P25989)
Applications	WB standard, antibody ELISA, Functional analysis, etc.
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Small delta antigen (Hepatitis delta virus/ISOLATE AMERICAN)

Note: