




H3(WR101)(Vaccinia Virus)

| | |
|----------------|--|
| Product Name | H3(WR101)(Vaccinia Virus) |
| Catalog No | IT-012-005Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea. |
| Protein | C-terminal 6xHis tagged H3(WR101)(Vaccinia Virus) Envelop protein (a.a. 1-284, GenBank# P07240) |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥90% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: H3(WR101)(Vaccinia Virus)

Note: