

# Anti-Spike RBD (CoVID-19/PAN)



Product Name	Anti- Spike RBD (CoVID-19/PAN)
Catalog No	IT-002-031-PAN-M1
Source	Mouse monoclonal antibody, clone 2H5G5
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Spike Protein from Novel Coronavirus (SARS-CoV-2/Wuhan/2019)(IT-002-032p).
Applications	ELISA, Neutralization, FACS, etc.
Epitope	Conformational
Specificity	Reacts with Spike S and RBD of multiple CoVID-19 variants as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

