

Anti-gp120/160(SIV/mac239)



Product Name	Anti-gp120/160(SIV/mac239)
Catalog No	IT-001-022
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	gp120(SIV/mac239) protein (Genebank accession #M33262)
Applications	Western Blotting (1:200-1:1000), Elisa (1:1000-1:5000)
Specificity	Reacts with SIV gp120/160 proteins.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

