




gp140 (HIV-1/Clade C)

Product Name	gp140 (HIV-1/Clade C)
Catalog No	IT-001-0031p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS only.
Protein	6xHis tagged HIV-1 gp140 (Clade C) protein (a.a.24-680) (Genebank #AY463234). Majority of protein exists as trimer.
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HIV-1 gp140 (Clade C) protein

Note: