

# Biotinylated gp120(Du151) (HIV-1/Clade C)



Product Name	Biotinylated gp120(Du151)(HIV-1/Clade C)
Catalog No	IT-001-139p-Biotin
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of biotinylated protein (0.5 mg/ml) in PBS+40% glycerol
Protein	6xHis tagged HIV-1 gp120 (Du151)(Clade C) protein (a.a.37-490) (Genebank No.AAL05313)
Applications	antibody ELISA, protein labeling, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

