

Biotinylated gp140 (SIV/mac239)



Product Name	Biotinylated gp140 (SIV/mac239)
Catalog No	IT-001-140p-Biotin
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of biotinylated protein (0.5 mg/ml) in PBS+40% glycerol
Protein	6xHis tagged HIV-1 gp140 (SIV/mac239) protein (a.a.20-692) (Genebank Accession No. M33262)
Applications	antibody ELISA, protein labeling, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

