

gp140(HIV-1/Bal.01)(Clade B)



Product Name	gp140(HIV-1/Bal.01)(Clade B)
Catalog No	IT-001-173ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gp140(HIV-1/Bal.01)(Clade B) (a.a.30-677) (Genebank Accession No. ABC55874)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

