

gp120 (C2-94UG114)(HIV-1/Clade D)



Product Name	gp120 (C2-94UG114)(HIV-1/Clade D)
Catalog No	IT-001-191p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged gp120 (C2-94UG114)(HIV-1/Clade D) (aa 1-496)(GenBank# AY494966)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

