


HA(Δ TM) (H7N9)(A/Shanghai/1/2013)



Product Name	HA(Δ TM) (H7N9)(A/Shanghai/1/2013)
Catalog No	IT-003-0073 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(Δ TM)(H7N9)(A/Shanghai/1/2013) (a.a.19-524) (GISAID #: EPI439486). This is the HA protein from the outbreak of the novel H7N9 virus in China since late March, 2013.
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H7N9)(A/Shanghai/1/2013)