



HA(Δ TM)(B/Brisbane/60/2008)

Product Name	HA(Δ TM)(B/Brisbane/60/2008)
Catalog No	IT-003-B3 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(B/Brisbane/60/2008) protein (a.a.16-548) (Genebank No. ACN29380)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required
Use Limitation	For research use only, not for use in diagnostic procedures.

