




Anti-GP (Lassa Virus)

Product Name	Anti-GP(Lassa Virus)
Catalog No	IT-008-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl rabbit anti-sera with 0.1% sodium Azide.
Immunogen	Glycoprotein protein of Lassa virus (Genebank No. AAL13212)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Glycoprotein protein of Lassa virus.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing GP(Lassa) protein

Note: Both cleaved and non-cleaved GP shown