


## matrix (Oz virus)



---

Product Name	matrix (Oz virus)
Catalog No	IT-025-009Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged matrix (Oz virus) (a.a. 1-270, GenBank# YP_009553283)
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: matrix (Oz virus)

Note: