

Anti- nucleocapsid protein [Mumps virus strain Jeryl Lynn]



Product Name	Anti- nucleocapsid protein [Mumps virus strain Jeryl Lynn]
Catalog No	IT-030-001M2
Source	Mouse monoclonal antibody, clone 3B10F2
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing nucleocapsid protein [Mumps virus strain Jeryl Lynn] (IT-030-001Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with nucleocapsid protein [Mumps virus strain Jeryl Lynn] as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

