

# Anti-M (hMPV)



Product Name	Anti-N (hMPV)
Catalog No	IT-031-002M1
Source	Mouse monoclonal antibody, clone 5B3E6
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant M(hMPV) (IT-031-002Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with M (hMPV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

