

# Human CD79b



Product Name	Human CD79b
Catalog No	IT-300-199p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human CD79b, namely B-cell antigen receptor complex-associated protein beta chain / Ig-beta / B-cell-specific glycoprotein B29 (aa 29-159; GenBank# P40259) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human CD79b

Note: