



NS1 (Usutu virus)

Product Name	NS1 (Usutu virus)
Catalog No	IT-006-0066p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100µg of purified protein (1mg/ml) in PBS.
Protein	C-terminal 6xHis tagged NS1(Usutu virus)(USUV) antigen (a.a. 794-1145)(Genbank Accession No. Q5WPU5).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

