

Human IBP1/IGFBP1/ Insulin-like growth factor-binding protein 1



Product Name	Human IBP1/IGFBP1/Insulin-like growth factor-binding protein 1
Catalog No	IT-010-035p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human IBP1/IGFBP1/Insulin-like growth factor-binding protein 1 (a.a. 26 - 259) (GenBank Accession No. P08833) with 6xHis tag at C-terminal.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95%, based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

