

# Env (SRV-8)



---

Product Name	Env (SRV-8)
Catalog No	IT-019-008p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Envelope protein (Simian retrovirus 8/SRV-8/SUZ/2012) (a.a. 23-516) (GenBank Accession No. AOS48101) with 6xHis tag at C-terminal. Majority are trimers.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95%, based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

