

glycoprotein gE (31-350)(HHV-3/VZV/Dumas)



Product Name	glycoprotein gE (31-350)(HHV-3/VZV/Dumas)
Catalog No	IT-005-043p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Varicella-zoster virus glycoprotein gE (31-350)(HHV-3/VZV/Dumas) (aa 31-350; GenBank# P09259) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: glycoprotein gE (31-350)(HHV-3/VZV/Dumas)

Note: