


Spike RBD (CoVID-19/BA.3.2 Cicada variant), His-tagged



Product Name	Spike RBD (CoVID-19/BA.3.2 Cicada variant), His-tagged
Catalog No	IT-002-036p-Cicada
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	The new coronavirus Spike protein RBD (CoVID-19/BA.3.2 Cicada variant) with 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: RBD (CoVID-19/BA.3.2), his-tagged

Note: